


# CHEMISTRY

## -LABORATORY INTRODUCTION-



---

## TABLE OF CONTENTS

---

LABORATORY INTRODUCTION.....	0
GENERAL RULES TO FOLLOW DURING AN EXPERIMENT .....	2
SAFETY WORKSHEET.....	4
LIST OF LABORATORY EQUIPMENT .....	7
GENERAL LABORATORY PROCEDURES WORKSHEET .....	9
HOW TO MEASURE LENGTH.....	10
HOW TO POUR A LIQUID .....	11
HOW TO MEASURE VOLUME .....	12
HOW TO LIGHT A BURNER.....	13

---

**LAB GENERAL RULES TO FOLLOW DURING AN EXPERIMENT<sup>1</sup>**

---

**1. KNOW WHAT YOU ARE DOING!**

- Read all directions before you start to work. Take notes.
- When in doubt, ask the teacher to explain.
- Use only the amounts and materials listed in the experiment.
- Make substitutes only when told to do so by your teacher.

**2. OBSERVE LABORATORY SAFETY PROCEDURES!**

- Always wear safety goggles and a laboratory apron while working in the laboratory.
- Fasten long hair back.
- Never taste, eat, or drink anything in the laboratory.
- Immediately wash off any chemical spilled on your skin. Flood the area with copious amounts of running water for at least 2 minutes. Tell, or have someone tell, the teacher immediately.
- If your skin itches or burns, flood the area of skin with water. Tell, or have someone tell, the teacher immediately.
- Keep your face back away from containers being heated, and don't lean over your work area.
- Avoid inhaling toxic fumes by using the ventilation hood.
- Wash your hands before leaving the laboratory.
- Report all accidents to your teacher at once.

**3. PREVENT ACCIDENTS FROM OCCURRING!**

- Work in a quiet, businesslike manner.
- Avoid moving about the laboratory.
- Clean up spilled chemicals by first flooding them with water and then wiping them up with a sponge.
- Keep the drawers and cupboard doors of your laboratory desk closed unless you are removing or putting something away.
- Don't reach over a burner. (Sometimes you can't see a flame!)
- Turn off a burner as soon as you are finished heating something, and remember the iron support ring stays hot for a while after the burner has been turned off.
- Put burned matches, broken glassware, and used chemicals in the places designated by your teacher.
- Know where the fire extinguisher, nearest fire alarm box, fire blanket, and eyewash device are located and how to use them.
- Report all accidents to your teacher immediately.

---

<sup>1</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York.

**4. AVOID CONTAMINATION OF CHEMICALS!**

- Never return any excess chemical to its stock bottle (original container). Discard the extra chemical, as you would used chemicals.
- Replace caps, stoppers and corks immediately after removing a chemical from its bottle. The chemical may be affected by the air or moisture.
- Use only clean, dry spoons, spatulas, or scoopulas to remove solids from their containers. Rinse and dry stirring rods before using them to stir other liquids.
- Be sure all glassware is clean before using. Drops of water will not cling to clean surfaces of glass.

**5. STUDY OBSERVATIONS BEFORE ANSWERING QUESTIONS!**

- Look for differences and similarities. Hunt for "patterns" in the data.
- Think what the question means before answering it. See if the textbook (or other references) gives you any help.
- Look up the meaning of any term you do not understand.
- Record any other questions or thoughts you may have that might lead to future experimentation or discussion.



"Always remember to write down your laboratory observations and data. You never know when you may need it."

LAB

**SAFETY WORKSHEET<sup>2</sup>****INTRODUCTION**

Because you use many chemicals in the chemistry laboratory, it is important to know something about them. Most chemicals are dangerous - e.g., poisonous when they enter the body, burn or irritate the skin, eyes and/or membranes of your digestive and respiratory tracts. Some chemicals destroy fabrics, especially synthetic or protein fibers. Chemicals can even make some things burn. Generally, most chemicals combine at a reasonable rate; however, under certain conditions they can get out of control. Knowing these things about chemicals helps you to prevent accidents and to work safely in the laboratory.

**QUESTIONS**

1. While you are working in the laboratory, why should you:
  - a. wear safety goggles? \_\_\_\_\_
  - b. fasten long hair back? \_\_\_\_\_
  - c. wear a laboratory apron? \_\_\_\_\_
  - d. use only the amounts and materials given in the directions?  
\_\_\_\_\_
  - e. use only clean glassware? \_\_\_\_\_
2. Why are unused chemicals never returned to their original containers?  
\_\_\_\_\_
3. Rings and other jewelry should not be worn in the lab. Why?  
\_\_\_\_\_
4. Why are bottles of chemicals kept tightly closed except when a quantity of chemical is being removed?  
\_\_\_\_\_
5. Why should the teacher be informed if you are wearing contact lenses?  
\_\_\_\_\_
6. Why are most chemicals stored in glass or plastic containers rather than metallic ones?  
\_\_\_\_\_
7. How should you dispose of chemicals that you used in the experiment?  
\_\_\_\_\_

<sup>2</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York. p. 6-8.

8. If your skin itches, burns, or appears red while you are using chemicals, what should you do? (a) \_\_\_\_\_  
(b) \_\_\_\_\_
9. What procedures should you follow to:
- report an accident resulting in injury to you or someone else? \_\_\_\_\_
  - report spilling a chemical or breaking equipment? \_\_\_\_\_
  - cope with a fire on your laboratory bench? \_\_\_\_\_
  - dispose of cracked or broken glassware? \_\_\_\_\_
10. Where are the following items located in the classroom?
- fire extinguisher: \_\_\_\_\_
  - fire blanket: \_\_\_\_\_
  - shower: \_\_\_\_\_
  - eye wash device: \_\_\_\_\_
11. Why are desk drawers and cupboard doors kept closed unless you are taking something out or putting it away?  
\_\_\_\_\_
12. Why is it important to work in a business-like manner, quietly, and with a minimum of distraction or migration around the classroom?  
\_\_\_\_\_
13. What is the meaning of the following words often found on the labels of chemicals?
- caustic: \_\_\_\_\_
  - corrosive: \_\_\_\_\_
  - hazardous: \_\_\_\_\_
  - volatile: \_\_\_\_\_
14. Why are you not allowed to eat or drink anything in the chemistry classroom?  
\_\_\_\_\_
15. Why should you read the label on a bottle twice before using its contents?  
\_\_\_\_\_
16. You should hold a bottle of liquid with the label toward the palm of your hand when you pour out the liquid. Why?  
\_\_\_\_\_
17. A cap, cork, or stopper from a bottle of liquid is not put down on a desk while you are pouring out some of the liquid. Why?  
\_\_\_\_\_

18. Explain the statement: "*You are only as safe as the LEAST safe person in the laboratory.*"















---

---

---

















# LAB LIST OF LABORATORY EQUIPMENT <sup>3</sup>

## LIST OF APPARATUS FOR STUDENT USE

DESCRIPTION	APPARATUS	USE	DESCRIPTION	APPARATUS	USE
glass common sizes 100 mL 250 mL 400 mL marked on the beaker	 400 mL <b>BEAKER</b>	as a container, like a cup may be heated	10 centimeter (cm) ruler, plastic divided into centimeter and millimeter (mm) divisions	 10 CM RULER	to measure length
glass common sizes 125 mL 250 mL 500 mL marked on the flask	 <b>ERLENMEYER FLASK</b>	may be heated	triangular wire frame with clay material coverings	 <b>PIPESTEM TRIANGLE</b>	to support the crucible
	 <b>FLORENCE FLASK</b>		small porcelain dish with cover	 <b>CRUCIBLE AND COVER</b>	to heat small amounts of solid material at high temperature
glass marked with a milliliter (mL) scale size divisions	 <b>GRADUATED CYLINDER</b>	to measure volume	hardened asbestos	 <b>ASBESTOS SQUARE</b>	to place under hot apparatus
50 mL 1.0 mL 35 mL 0.2 or 10 mL 0.1 mL			wire screen asbestos center	 <b>WIRE GAUZE</b>	to spread the heat of a flame
glass several sizes	 <b>TEST TUBE</b>	many uses can be heated	metal heating device connected to gas outlet with rubber tubing	 <b>LABORATORY BURNER</b>	to heat chemicals in beakers or test tubes
metal clamp with a spring handle	 <b>TEST TUBE CLAMP</b>	to hold a test tube	metal rod upright heavy base	 <b>RING STAND</b>	a support with many uses
metal	 <b>TONGS</b>	to pick up and hold apparatus			

<sup>3</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York. p. 4-5.



DESCRIPTION	APPARATUS	USE	DESCRIPTION	APPARATUS	USE
iron ring with screw fastener several sizes	 <b>IRON RING</b>	to fasten to the ring stand as a support for apparatus	brush with wire handle	 <b>TEST TUBE BRUSH</b>	to scrub glass apparatus
metal clamp with 1. screw fastener 2. swivel and lock nut 3. adjusting screw 4. curved clamp	 <b>BURET CLAMP</b>	to hold apparatus may be fastened to the ring stand	glass rod	 <b>STIRRING ROD</b>	to stir combinations of materials to use in pouring liquids
heavy porcelain dish with grinder	 <b>MORTAR AND PESTLE</b>	to grind chemicals to a powder	porcelain dish	 <b>EVAPORATING DISH</b>	as a container for small amounts of liquid being evaporated
may be of metal or porcelain	 <b>SPATULA</b>	to transfer solid chemicals in weighing	thick glass	 <b>GLASS PLATE</b>	many uses (should not be heated)
metal file with three cutting edges	 <b>TRIANGULAR FILE</b>	to scratch glass to file	curved glass	 <b>WATCH GLASS</b>	may be used as a beaker cover may be used in evaporating very small amounts of liquid
short length of rubber tubing	 <b>RUBBER CONNECTOR</b>	to connect parts of apparatus	glass or plastic	 <b>FUNNEL</b>	to hold a filter paper may be used in pouring
metal clamp with finger grips	 <b>PINCH CLAMP</b>	to clamp a rubber connector	glass tip with rubber bulb	 <b>MEDICINE DROPPER</b>	to transfer small amounts of liquid
rack; may be wood, metal or plastic	 <b>TEST TUBE RACK</b>	to hold test tubes in an upright position	metal	 <b>FORCEPS</b>	to pick up or hold small objects

---













**LAB GENERAL LABORATORY PROCEDURES WORKSHEET<sup>4</sup>**


---

**QUESTIONS**

1. Name the apparatus or glassware is used?
  - a. to pick up and hold apparatus? \_\_\_\_\_
  - b. to clean glassware? \_\_\_\_\_
  - c. as a cover for a beaker? \_\_\_\_\_
  - d. to transfer solid chemicals in massing? \_\_\_\_\_
  - e. as a container for small amounts of liquid being evaporated? \_\_\_\_\_  
or \_\_\_\_\_
  - f. to grind chemicals to a powder? \_\_\_\_\_
  - g. to measure a volume of a liquid? \_\_\_\_\_
  - h. to measure length? \_\_\_\_\_
  - i. to measure mass? \_\_\_\_\_

2. What is the name for each of the following pieces of laboratory equipment?  
(N.B. Objects are not shown to scale.)



a. 	b. 	c. 
d. 	e. 	f. 
g. 	h. 	i. 
j. 	k. 	l. 

<sup>4</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York. p. 8-10.

---

**LAB HOW TO MEASURE LENGTH<sup>5</sup>**


---



FINAL READING ( $R_2$ ) MINUS INITIAL READING ( $R_1$ ) = LENGTH

$$R_1 = 1.0 \text{ cm} \quad R_2 = 7.4 \text{ cm}$$

$$7.4 \text{ cm} - 1.0 \text{ cm} = 6.4 \text{ cm LONG}$$

---

<sup>5</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York. p. 23


---

**LAB HOW TO POUR A LIQUID<sup>6</sup>**


---

Remove the cap to the bottle. Do not put cap on the counter in a way that would allow liquid on the counter top – when you do this, you have essentially spilled chemical onto the countertop. Moreover, when you replace the cap on the bottle, you can contaminate the liquid if there was anything on the countertop.

Shown below is the correct way to remove the stopper to a bottle containing a liquid. A similar procedure should be used if the cap is a screw-cap, parafilm<sup>®</sup> wrapping, etc.



Generally, liquids can be poured carefully directly into the desired container. However, one common technique, used especially for pouring liquids that also contain suspended solids, is **DECANTING**, which is shown below.



Also, don't hesitate to use a funnel when pouring into a small-mouth container.




---

<sup>6</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York. p. 17

## LAB HOW TO MEASURE VOLUME<sup>7</sup>

1. Decide what each unmarked line represents:



A. On the above scale (Figure 1):


1)

$$\text{Divide: } \frac{\text{difference}}{\text{number of lines}} = \frac{\text{amount}}{1 \text{ line}}$$

2)



$$\frac{5 \text{ mL}}{5 \text{ lines}} = 1 \text{ mL per line}$$

B. Look carefully at the scale. Some graduated cylinders are marked differently.



C. Read at the bottom of the meniscus (the curved surface of the liquid).

Hold the graduated cylinder at eye level so that it is high enough that you can't see the back of the nearest number.



<sup>7</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York. p. 22

---

**LAB HOW TO LIGHT A BURNER<sup>8</sup>**


---

**OVERVIEW:**

Fire requires three parts: (1) heat, (2) fuel, and (3) oxygen. Remove any one of these and the fire will be extinguished. There are several styles of burners available but all share the combining of heat, fuel, and oxygen. Gas intake can be adjusted at the main gas valve or by a valve at the base of the burner. Oxygen intake can be adjusted at the air intake. The style of the burner and the source of heat will determine the steps you follow. Listed below are the two most common burners and two ways for lighting either type.

The hottest part of any flame is the tip of the inner blue cone.




Figure 1.

**Fischer Burner**

1. Use rubber tubing to connect the burner to the gas valve. Make sure that the chimney is lowered and just snug (no he-man tightening: it will strip the threads), or that the collar is closed. This will make it easier to light the gas.
2. Turn on the main gas valve.
3. Strike the striker and ignite the gas. (If you smell gas, turn off the main gas valve and wait a couple of minutes for the odor to dissipate.)
4. Adjust the height of the flame with the burner gas valve (Figure 2) – again, just so it is snug or you will damage the valve.
5. Adjust the air intake to obtain a blue cone (Figure 1).




Figure 2.

**Bunsen Burner**

1. Use rubber tubing to connect the burner to the gas valve. Make sure that the chimney is lowered and snugly tight (no he-man tightening: it will strip the threads), or that the collar is closed. This will make it easier to light the gas.
2. Strike match.
3. Turn on the main gas valve.
4. Bring lighted match up to the chimney to ignite the gas. (If you smell gas, turn off the main gas valve and wait a couple of minutes for the odor to dissipate.)
5. Adjust the height of the flame with the main gas valve. Adjust the air intake (by raising/lowering chimney or by adjusting the air intake) to obtain a blue cone (Figure 1).




Figure 3.

---

**NOTES:**


---

1. **Remember: hot labware looks the same as cold labware!**

---

<sup>8</sup> Bolton, Ruth P., Lamphere, Elizabeth V., and Menesini, Mario. (1979). *Laboratory Experiments in ACTION CHEMISTRY*. Holt, Rinehart and Winston, Publishers, New York. p. 22